

Intoxilyzer - Alcohol Analyzer
Model 9000

SN: 90-003795

Test Date: 06/13/2025

Test Time: 04:47

Location: -----

10701 S ESTRELLA PKY

Software Version: 9439.01.00

Standard Lot #: 2842100A1

Standard Value: 0.100

Expiration Date: 11/05/2026

Bottle #: 24

Last Changed By: MIYAZATO, KORY

Permit #: 30215

Operator: MIYAZATO, KORY

Permit #: 30215

Agency: GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	04:47:26
Calibration Chk	0.100	04:47:48
Air Blank	0.000	04:48:28

No RFI Detected

Successfully Completed Test Sequence

Intoxilyzer - Alcohol Analyzer

Model 9000

SN: 90-004605

Test Date: 06/13/2025

Test Time: 04:28

Location:

Software Version: 9439.01.00

Standard Lot #: 302-402972207

Standard Value: 0.100

Expiration Date: 02/20/2026

Bottle #: 19

Last Changed By: MIYAZATO, KORY

Permit #: 30215

Operator: MIYAZATO, KORY

Permit #: 30215

Agency: GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	04:28:32
Calibration Chk	0.098	04:28:53
Air Blank	0.000	04:29:33

No RFI Detected

Successfully Completed Test Sequence

EXHIBIT I-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
INTOXILYZER MODEL 9000

PERIODIC MAINTENANCE AND STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST K Miyazato / 1139 AGENCY GOODYEAR
DATE 6/13/25 TIME 0500
INTOXILYZER SERIAL # 90-003795

- ☒ 1. Ensure that gas tank is attached and contains a standard alcohol concentration -100 AC

DIAGNOSTIC TESTS

- ☒ 1. Clock time check
☒ 2. Date check

OPERATIONAL TESTS

- ☒ 1. Deficient Subject Test (Proper Sample Recognition)
Deficient Sample printed
☒ 2. Alcohol-free Subject Test (Proper Sample Recognition):
0.000 AC
☒ 3. Mouth Alcohol Subject Test (Proper Sample Recognition):
Invalid Sample – Begin new deprivation period printed
☒ 4. Radio Frequency Interference Test (Error Recognition):
RFI Detect printed
☒ 5. Standard Calibration Check:
0.100 AC
☒ 6. Air Blanks Completed
☒ 7. Timer Reset

Not a Successfully Completed Test Sequence will be printed

Instrument is operating properly and accurately. YES ☒ NO ☐

COMMENTS:

SIGNATURE K Miyazato / 1139

ARIZONA FORENSIC BREATH ALCOHOL ANALYTICAL REPORT

11 N 145TH AVE

INSTRUMENT INFORMATION

Analytical Instrument: Intoxilyzer 9000
Serial Number: 90-003795
Software Version: 9439.01.00
Analytical Report Number: 37950613250720

QAS: MIYAZATO, KORY
QAS Permit #: 30215
Agency: GOODYEAR PD
Last 31-Day Check: 05/14/2025
Last Annual Maintenance: 10/02/2024

SUBJECT INFORMATION

Name: 31-Day Check
Test Date: 06/13/2025
Driver's License #:
Gender:
Date of Birth:
Age:

Weight:
State of Issue:
Driver's License Expiration:
Deprivation Start Time: 04:40
15 - Minute Deprivation: Yes

OPERATOR INFORMATION

Name: MIYAZATO, KORY
Agency: GOODYEAR PD

Permit #: 30215

STANDARD INFORMATION

Standard Value: 0.100
Standard Lot #: 2842100A1
Expiration Date: 11/05/2026
Bottle #: 24
Last Changed By: MIYAZATO, KORY

Permit #: 30215

Test	g/210L	Time
Air Blank	0.000	05:00:41
Diagnostic Test	PASS	05:01:16
Air Blank	0.000	05:01:54
Calibration Chk	0.100	05:02:15
Air Blank	0.000	05:02:55
Subject Test 1	DEF*	05:06:15
Air Blank	0.000	05:06:55
Air Blank	0.000	05:07:35
Subject Test 2	0.000	05:07:55
Air Blank	0.000	05:08:36
Wait		05:12:40
Air Blank	0.000	05:13:17
Subject Test 3	INV**	05:13:35
Air Blank	RFI***	05:13:48

*Deficient Sample
**Invalid Sample - Begin new deprivation period
***RFI Detect

RESULTS

EXCEPTION MESSAGES

* Deficient Sample
** Invalid Sample - Begin new deprivation period
*** RFI Detect

OPERATOR COMMENTS

Not a Successfully Completed Test Sequence

All air blank results must be 0.000.
Consecutive subject test results must not differ by more than 0.020 g/210L.
Standard check results must be $\pm 10\%$.

Intoxilyzer - Alcohol Analyzer
Model 9000

SN: 90-003795

Test Date: 06/13/2025

Test Time: 05:14

Location: 11 N 145TH AVE

Software Version: 9439.01.00

Standard Lot #: 2842100A1

Standard Value: 0.100

Expiration Date: 11/05/2026

Bottle #: 24

Last Changed By: MIYAZATO, KORY

Permit #: 30215

Operator: MIYAZATO, KORY

Permit #: 30215

Agency: GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	05:14:39
Calibration Chk	RFI*	05:14:58
Air Blank	RFI**	05:15:19

*RFI Detect

**RFI Detect

Intoxilyzer - Alcohol Analyzer
Model 9000

SN: 90-003795

Test Date: 06/13/2025

Test Time: 05:15

Location: 11 N 145TH AVE

Software Version: 9439.01.00

Standard Lot #: 2842100A1

Standard Value: 0.100

Expiration Date: 11/05/2026

Bottle #: 24

Last Changed By: MIYAZATO, KORY

Permit #: 30215

Operator: MIYAZATO, KORY

Permit #: 30215

Agency: GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	05:15:58
Calibration Chk	0.100	05:16:20
Air Blank	0.000	05:17:00

No RFI Detected